

Installing AORTIC upgrade – from version 7.0.0 (or higher)

IMPORTANT

In order to use the new Admission and Discharge Summary (ICU Adms Form) you must have the AORTIC program files installed in the following directory: C:\Program Files\AORTIC.

Download the latest version 9.2.1 upgrade software from the ANZICS website or obtain a copy of the upgrade for CD.

If you are currently using any version of AORTIC 7.x.x or higher follow the upgrade instructions below:

Before you begin installation, perform a backup of the AORTIC Database and exit from every open application so that only Windows is running. There should be no applications listed in the task bar or in floating tool bars.

For a network environment, the update must be run on each client copy of AORTIC. All end users need to exit the program before commencement of the upgrade. All client applications need to be upgraded before data entry is recommenced.

Note: Installation of the update replaces the front end program files in an existing AORTIC installation. The back end database files (installed on a standalone PC or network) containing patient information are not reloaded.

To install the AORTIC database program upgrade:

1. Insert the update AORTIC CD into your CD-ROM Drive. Select the .exe file. The main WinZip dialog box will appear (Figure 1.13)




Figure 1.13 WinZip® main dialog box

2. The default location for installation of the files is C:\Program Files\AORTIC. Click **Unzip** to install the files here.

3. If the files are to be located other than C:\Program Files\AORTIC, click on **Browse**. A dialog box to enable selection of the required folder will be displayed (Figure 1.14). Folders, or directories and subdirectories can be seen by clicking on the + symbols. When the chosen folder is selected click the **OK** button. The main WinZip dialog box will appear again. Click **Unzip**.




Figure 1.14 Browse Folder

4. An indicator bar will appear at the bottom of the dialog box to indicate progress. A message box (Figure 1.15) will appear to indicate the successful installation of the files. Click **OK** then Close.



Figure 1.15 WinZip Message showing success